



Aged Care



Aged Care

We design and manufacture a range of specialised patient care equipment for the aged care sector, working with both head office and at branch levels to ensure best outcomes.



I-move



I-move Drawer



I-move Slim



Contour Recline



Contour Energy



Aquarius



Aquatuff



Aquatuff Bari



Aquatuff Electric



Change Table



Custom Liner Colours

I-move

A height adjustable rounds trolley with storage that is light and easy to move in tight spaces.



Key Benefits

Compact & simple: Light, easy to move and unobtrusive in tight spaces is the perfect partner to make sure your nursing rounds are efficient and easy.

Tech support: Upgrade your I-move to support any and all of your IT requirements required for the job.

Safe & secure: Option to add a Kensington lock to restrain valuable laptop hardware to keep safe as well as cable management options.

Why Buy?

The i-move is an affordable and height adjustable rounds trolley. When it comes to moving laptops, tablets or simply pens and paper along with other clinical items, the only practical solution has been the traditional metal dressing trolley.

How good would it be to have a trolley that is light and easy to move, unobtrusive in tight spaces and that you can comfortably sit or stand at? Coupled with this is a convenient storage shelf and our ability to customise the top size and shape to totally suit your requirements!

I-move



Storage space

Convenient storage space with retaining lip across the front.



Customisable

Use our standard work surface or customise to suit your specific requirements.



Height adjustable

Simply squeeze the lever under the top and effortlessly adjust to your desired work height.



Storage compartment

Add a removable stainless steel storage compartment for cables / pens etc.



Power supply holder

Securely hold a laptop power supply for easy charging.



Kensington lock

Restrain valuable laptop hardware to keep it safe.



Add your IT

Add your own holders to suit the specific IT hardware in use.

Specifications

Size: Compact, easy to clean top - 550mm wide x 560mm deep

Height adjustable: From 840mm to 1140mm

Storage size: Convenient 100mm storage space with lip

Product Codes

I-Move Trolley; MCT-IM-LHA

Removable Storage Container - MCT-IMA-SC

Power Supply Holder - MCT-IMA-PCH

Kensington Lock - MCT-IMA-KL

I-move Drawer

A height adjustable rounds trolley with a secure locking drawer.



Key Benefits

Safe & secure: Quality metal lockable drawer with ball bearing runners allows you to safely store valuables and secure them.

Custom design: Choose between a standard work surface or customise to suit your desired requirements.

Tech savvy: Upgrade your I-move Drawer to support a wide range of IT hardware.

Why Buy?

This i-move is a height adjustable rounds trolley with a secure locking drawer. When it comes to moving laptops, tablets or simply pens and paper along with other clinical items, the only practical solution has been the traditional metal dressing trolley.

How good would it be to have a trolley that is light and easy to move, unobtrusive in tight spaces and that you can comfortably sit or stand at? Coupled with this is a convenient storage drawer and our ability to customise the top size and shape to totally suit your requirements!

I-move Drawer



Lockable drawer

Safe 400 wide x 390 x 150mm deep drawer on quality ball bearing runners.



Compact design

Move your items easily and safely with this practical rounds trolley.



Customisable

Use our standard work surface or customise to suit your specific requirements.



Height adjustable

Simply squeeze the lever under the top and effortlessly adjust to your desired work height.



Accessories for IT

Add holders or locks for a wide range of IT hardware.

Specifications

Size: Compact, easy to clean top - 550mm wide x 645mm deep

Drawer Size: 400 wide x 390 x 150mm deep drawer size

Height Range: Height adjustable from 930mm to 1230mm

Options: Laptop security lock and tablet mount systems available

Product Code

MCT-IM-LHA-LD

I-move Slim

No under top storage for an ultra lightweight, height adjustable and easy to move work surface.



Key Benefits

Compact & simple: Lightweight design creates a practical rounds trolley for people 'on the go'

IT savvy: Add holders or locks for a wide range of IT hardware. Go to the I-Move IT trolley to see all options.

Custom design: Choose between a standard work surface or customise to suit your desired

Why Buy?

How good would it be to have a height adjustable, compact work surface that is very light and easy to move in tight spaces and that you can sit or stand at? For any area of any healthcare facility, this trolley brings a very useful work surface with you, wherever you go.

I-move Slim



Ultra thin

A lightweight, practical rounds trolley for people 'on the go'.



Customisable

Use our standard work surface or customise to suit your specific requirements.



Height adjustable

Simply squeeze the lever under the top and effortlessly adjust to your desired work height.



Accessories for IT

Add holders or locks for a wide range of IT hardware. Go to the I-Move IT trolley to see all options.

Specifications

Size: Compact, easy to clean top - 550mm wide x 560mm deep

Height Range: Height adjustable from 730mm to 1030mm

Options: Laptop security lock system, tablet mount & rear cable storage compartment options

Product Code

MCT-IM-LHA-ST

Contour Recline

A versatile chair to streamline medical examination, procedures and recovery across many healthcare facilities.



Key Benefits

Save Space: Fit more chairs in the same space by increasing throughput with faster recoveries from patients sitting up right.

Hygiene conscious: Smooth, gap free platform aids cleaning and helps prevent fluid ingress to meet the highest infection control standards.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with quality, care and customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance and fast guaranteed parts.

Application

Ideal for Day Surgery, Recovery, Transit Lounge, Fast Track ED & Cardiac

Why Buy?

Make better use of your valuable space and streamline your department with Modsel's compact procedure chairs. They are easy to move in and out of tight areas meaning a higher patient throughput and reduced waiting times. Patients can be positioned easily from sitting upright to a full recline, or lay completely flat with the stability of twin columns under the chair which can be used for CPR or specialised procedures. We've taken care with every aspect of design to ensure comfort and efficiency for both the patient and carer. Pop-up push handles ensure best operator positioning while preventing transfer movement to the patient.

Contour Recline



Low and safe

Excellent 460mm low height with fold-under drop-sides.

The ultra-low height is combined with drop-sides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patients legs or clothes.



Easy egress

Forward tilt facilitates getting on and off.

For aged or disabled patients, a gentle incline and the long drop-side armrests may assist in safer transfer.



Lies flat

The backrest and leg rest move independently or together from chair to bed.

You now have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Smart transfer

Conveniently stored handles unfold for staff friendly transportation.

Lock in the double steer system, unfold the push handles and move forward effortlessly on the large castors. The unique handle design correctly positions the carer while preventing unwanted forces being applied directly to the patient.



Wide option

100mm wider top option available.

Change the seat width from 650mm to 750mm wide to easily accommodate larger patients. Use our range of colour coding choices if required to assist in quick identification.



Clear handset

Large simple functions with colour assistance.

Position the patient quickly and easily with the touch of a button - either by independently moving sections or applying one touch contouring.

Specifications

Size: 780 x 1150mm overall (1750mm lying flat)

Seat width: 650mm with 100mm wider option available with 100mm thick premium, removable mattress

Height Range: 460 - 840mm with lie-flat and forward / reverse tilt

SWL: 300kg

Castor Size: Twin wheel 150mm castors with double steer

Package Includes: Standard IV, vertical oxy holder, dropdown footrest & meal tray

Product Codes

Contour Recline: MPC-C-RET

Contour Recline Package MPC-C-RET-P

750mm wide top option: MPC-CA-WCT

Contour Energy

A transport chair with ZERO effort, 4-way power drive reducing the risk of strain during patient transport.



Key Benefits

Ergonomics & intuitive: Drive effortlessly in any direction with only the slightest of touch, paired with perfect positioning of the drive controls for ease of use.

Easy access: The large footplates offer generous support to your patients and lift quickly for clear egress.

Powerful movement: A rotating electric drive wheel provides excellent traction, while center steering makes patient transport easier over longer distances.

Why Buy?

Take all the human effort out of pushing and maneuvering chairs with an intuitive and intelligent drive system designed for easy use by any operator. Natural movement controls the speed and direction of the chair with the ability to stop very quickly by pulling back. Help prevent the most common carer strain injuries by taking the back twist out of moving heavy weight forwards, backwards and sideways. Patients are positioned comfortably with a gentle recline and feet fully supported. We've taken care with every aspect of design to ensure comfort and efficiency for both the patient and carer. Customise for a wide range of healthcare areas with relevant accessories.

Contour Energy



Ergonomic

Drive effortlessly in any direction by only light touch.



Intuitive

Perfect positioning of the drive controls for the operator.



Support

The large foot plates offer generous support to the patient.



Easy access

The large foot plates lift quickly and easily for clear egress.



Full clearance

Complete frontal clearance by lifting and rotating foot plates outwards.



Removable armrests

The side armrests can be fully removed for multi direction egress.



Versatile

Use the storage shelf during transport or customise with accessories such as IV poles, oxy holders and seat belt.

Specifications

Size: 750mm W x 1300mm L x 1180mm H

Seat Size: 650mm W x 520mm D x 550mm H

SWL: 300kg

Castor Size: 150mm

Package Includes: standard IV & vertical oxy holder

Product Code

Contour Energy: MPC-C-ENT-EDU

Contour Energy Package: MPC-C-ENT-EDU-P

Aquarius

A standard shower trolley designed to meet budget constraints.



Key Benefits

Safe & secure: The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfers.

Hygiene conscious: Tilt the top 6 degrees and lower the drain hose tip to quickly remove waste water.

Custom adjustments: Raise or lower the height of your trolley to suit bed transfer or best working positions. You can also lock the removable end rails into two positions, allowing 120mm extra space for longer patients.

Application

Designed for standard bathroom use.

Why Buy?

We are able to bring our Modsel shower trolley design quality to you at an affordable price, helping to meet tight budget constraints and maximising the number of working units available.

The refined framework and proven mechanics awards long life value and consists of considered, easy to clean surfaces coupled with quality twin wheel castors to ensure the highest level of maneuverability.

Aquarius



Safe transfer

The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfer.



Easy drainage

Tilt the top 6 degrees and lower the drain hose tip to quickly remove waste water.



Height adjustment

Raise or lower the height of your trolley to suit bed transfers or best working positions.



Length adjustment

Lock the removable end rails (on each end of the top) into two positions, allowing 120mm extra space for longer patients.



Junior size option

A special junior sized shower trolley developed for small patients. Working with children or small adults may not require a full sized shower trolley. Our junior sized version is also easier to move around in tight facilities and for storage.

Specifications

Size: 795 W x 1950mm including bumpers

Liner Inner Base Size: 550 W x 1700 L x 180mm Deep

Standard Hose Length: 1100mm, custom lengths available

SWL: 180kg

Height Range: 680 to 880mm

Castor Size: Quality 125mm twin wheel castors with total lock and steer

Construction: Powder coated steel

Junior Option: 795 W x 1680mm including bumpers

Junior Liner Base Size: 550 W x 1450 L x 180mm Deep

Product Codes

Aquarius - MST-AQM

Aquarius Junior - MST-AQM-J

Aquatuff

A stainless steel shower trolley with hydraulically operated 300kg SWL and wide top option.



Key Benefits

Safe & secure: The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfers.

Hygiene conscious: Tilt the top 2 degrees and lower the drain hose tip to quickly remove waste water.

Custom adjustments: Raise or lower the height of your trolley to suit bed transfer or best working positions. You can also lock the removable end rails into two positions, allowing 120mm extra space for longer patients.

Application

Designed for heavy use and demanding environments.

Why Buy?

Having equipment that is robust and completely reliable is of the utmost importance when caring for larger or heavier patients. We have over-engineered the componentry strength in this design and our decision to use total stainless steel reinforces the lifetime of value that this bariatric shower trolley will give. The hydraulic operation helps to safely lift 300kg without overload, whilst four twin wheels with brakes and wider top all increase versatility in special needs.

Aquatuff



Safe transfer

The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfer.



Easy drainage

Tilt the top 2 degrees and lower the drain hose tip to quickly remove waste water.



Height adjustment

Raise or lower the height of your trolley to suit bed transfers or best working positions.



Length adjustment

Lock the removable end rails (on each end of the top) into two positions, allowing 120mm extra space for longer patients.



Wide top option

Add a 100mm wider top for larger patients with wide shoulders.

Specifications

Size: 795 W x 1950mm including bumpers

Liner Inner Base Size: 550 W x 1700 L x 180mm Deep

Standard Hose Length: 1100mm, custom lengths available

SWL: 300kg (electric model available)

Height Range: 660 to 900mm

Castor Size: Quality 125mm twin wheel castors with total lock and steer

Construction: Powder coated stainless steel

Junior Option: 795 W x 1680mm including bumpers

Junior Liner Inner Size: 550 W x 1450 L x 180mm Deep

Wider Comfort Option: 100mm (internal min width of 660mm)

Product Code

Aquatuff Hydraulic; MST-AQT

Aquatuff Hydraulic Junior; MST-AQT-J

100mm wider top option - MST-SA-EWT

Aquatuff Bari

A robust, fully electric 300kg SWL shower trolley with extra wide base and top for Bariatric patients.



Key Benefits

Safe & secure: The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfers.

Hygiene conscious: Tilt the top 2 degrees and lower the drain hose tip to quickly remove waste water.

Custom adjustments: Raise your patient to a comfort position with a 15 degree electric backrest option available.

Application

Designed for heavy use and demanding environments.

Why Buy?

Having equipment that is robust and completely reliable is of the utmost importance when caring for larger or heavier patients. We have over-engineered the componentry strength in this design and our decision to use total stainless steel reinforces the lifetime of value that this bariatric shower trolley will give. The electric operation helps to safely lift 300kg without overload, whilst four twin wheels with brakes with wider base and top all increasing versatility in special needs.

Aquatuff Bari



Safe transfer

The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfer.



Easy drainage

Tilt the top 2 degrees and lower the drain hose tip to quickly remove waste water.



Backrest option

Raise your patient to a comfort position with a 15 degree electric backrest.



Length adjustment

Lock the removable end rails (on each end of the top) into two positions, allowing 120mm extra space for longer patients.

Specifications

Size: 1195 W x 1950mm including bumpers

Liner Inner (min) Base Size: 950 W x 1700 L x 180mm Deep

Standard Hose Length: 1100mm, custom lengths available

SWL: 300kg

Height Range: 680 to 900mm

Castor Size: Quality 125mm twin wheel castors with total lock and steer

Electric Backrest Option: 15 degree

Construction: Powder coated stainless steel

Product Code

Aquatuff Bari; MST-BARI

Electric backrest option - MST-SA-EBB

Aquatuff Electric

A stainless steel shower trolley with electrically operated 300kg SWL, wide top and backrest options.



Key Benefits

Safe & secure: The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfers.

Hygiene conscious: Tilt the top 2 degrees and lower the drain hose tip to quickly remove waste water.

Custom adjustments: Raise or lower the height of your trolley to suit bed transfer or best working positions. You can also lock the removable end rails into two positions, allowing 120mm extra space for longer patients

Application

Designed for heavy use and demanding environments.

Why Buy?

Having equipment that is robust and completely reliable is of the utmost importance when caring for larger or heavier patients. We have over-engineered the componentry strength in this design and our decision to use total stainless steel reinforces the lifetime of value that this bariatric shower trolley will give. The electric operation helps to safely lift 300kg without overload, whilst four twin wheels with brakes, electric backrest and wider top options all increase versatility in special needs.

Aquatuff Electric



Standard design

The standard top has a full length dropside with an electric backrest raise option for extra comfort if required.



Safe transfer

The dropsides fold underneath the top and the liner side lays down to provide gap free patient transfer.



Easy drainage

Tilt the top 2 degrees and lower the drain hose tip to quickly remove waste water.



Height adjustment

Raise or lower the height of your trolley to suit bed transfers or best working positions.



Length adjustment

Lock the removable end rails (on each end of the top) into two positions, allowing 120mm extra space for longer patients.



Wide top option

Add a 100mm wider top for larger patients with wide shoulders.

Specifications

Size: 795 W x 1950mm including bumpers

Liner Inner Base Size: 550 W x 1700 L x 180mm Deep

Standard Hose Length: 1100mm

SWL: 300kg

Height Range: 680 to 900mm

Castor Size: Quality 125mm twin wheel castors

Construction: Powder coated stainless steel

Junior Size: 795 W x 1680mm

Junior Liner Inner Base Size: 560 W x 1450 L x 180mm Deep

Wider Comfort Option: 100mm

Product Code

Aquatuff Electric; MST-AQTE

Aquatuff Electric Junior; MST-AQTE-J

Electric backrest option (excludes Junior size) - MST-SA-EBA

100mm wider top option (excludes Junior size) - MST-SA-EWT

Change Table

A purpose built, non-slip change mat on a sturdy, height adjustable and safe patient platform.



Key Benefits

Safe & secure: The fold up dropsides and removable end rails ensure high levels of patient safety when required.

Customisable design: Tailor to your requirements with options such as electric high/low, 300kg SWL and wide tops.

Why Buy?

This change mat expands the use of the Aquarius trolley into a very safe, height adjustable change table with multi-position, removable end rails and full length dropsides that fold back underneath the trolley top for ample work clearance. The mat securely fastens underneath the trolley top to prevent sliding and movement. This trolley is designed for harsh use and able to be rigorously cleaned.

Change Table



Safe design

The fold up dropsides and removable end rails ensure a high level of patient safety when required.



Customise

Tailor to your requirements with options such as electric high/low, 300kg SWL and wide tops.

Specifications

Size: 795 W x 1950mm including bumpers and 25mm high density foam core

Patient Surface Size: 600 W x 1820mm Long

Height Range: 60 to 900mm (depending on trolley model)

SWL: Hydraulic 180kg or electric 300kg

Castor Size: Quality 125mm twin wheel castors with total lock and steer

Options: Both junior and wider sized available

Product Code

MST-SA-CTU

Custom Liner Colours

Customise the colour of your shower trolley liner from the standard Royal Blue.

Key Benefits

- Heavy duty washable, non-stick vinyl
- Ask for the available colours

Product Code

- MST-SA-CLC

Why Buy

There are a selection of light and dark colours available in the high grade vinyl used for the Aqua liners, that you can choose as an optional upgrade. Please contact us to discuss your requirements and view a current colour chart.



Royal Blue

Notes

Notes

Notes
