

# Knee Break

With electric or manual operated options, the knee section raises for patient comfort during recovery or long stay.

## Key Benefits

- Fits both standard and x-ray Contour stretchers
- Minimal imaging coverage loss with unrestricted cassette travel
- Excellent 150mm raise height for a range of comfort options

## Product Code

- Electric - MPT-CA-KBE
- Manual - MPT-CA-KBH

## Why Buy

Add this option to the Contour Portare and Portare-X or twin column stretchers such as the Contour Classic and Multi-X to enhance patient comfort and recovery. The unique knee break system also fits to emergency imaging models and maintains maximum cassette travel freedom and vision.



Electric Knee Break



Manual knee break



Fold away handle

v-14012022