



X-Ray Imaging



X-Ray Imaging

Reliable and practical surfaces for imaging is critical because it increases patient throughput and improves operating costs, helping to keep you in control. Whatever the procedure, our stretchers are carefully designed to maintain this balance with the best patient access and effective use of imaging equipment.



Contour Portare-X



Contour Multi-X



Contour Barituff-X



Contour Endosurge



Contour Endo-X

Contour Portare-X

An ultra-low electric trauma stretcher with totally unobstructed x-ray clearance.



Key Benefits

Unobstructed x-rays: Supported by just four outside corners, our unique x-ray system gives total freedom of movement for cassettes - even with an electric backrest!

Safe & low: The 430mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting.

Smooth & steady: Operated by a very light touch pedal action, this wheel positively pre-loads onto the floor to ensure strong grip and rides up and down effortlessly over bumps and hollows to keep traction.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for ED & Imaging

Why Buy?

Helping patients on and off a trauma stretcher is safer and easier when the bed can be lowered down to 430mm, around the standard office chair height! This stretcher can feature our unique InfineX* system which allows totally unobstructed head to toe clearance for effortless and efficient x-ray cassette positioning. Smoothly moving patients helps improve their experience and our E-Point steering delivers a strong grip to the floor with a solid, perfectly centered wheel to allow the stretcher to always accurately track and turn. The Modsel x-ray stretchers are designed to be a reliable, cost-effective performer, no matter what the shift brings.

*Patent pending.

Contour Portare-X



Unobstructed x-ray tunnel

100% x-ray clearance of the entire body. Supported by just four outside corners, our unique patented InfineX system gives total freedom of movement for x-ray cassettes - even with an electric backrest! There are two styles available, another that allows the backrest to raise separately to the lower cassette support surface.



Ultra-low height

Safer and easier for unwell patients to get on and off. The 430mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Excellent high height

Quickly raise a patient up to the best working height. Lifting the platform height up to 880mm helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Head and foot tilt

A generous 15 degrees top tilt both directions. Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



E-Point steering

A perfectly centred steer wheel ensures the stretcher always accurately tracks and turns. Operated by a very light touch pedal action, this wheel positively pre-loads onto the floor to ensure strong grip and rides up and down effortlessly over bumps and hollows to keep traction. Turning left or right from either end gives the same steering result.



Easy customisation

Increase top width, change colours, mattress and dropside styles. Every trauma department has their own unique requirements to be able to function safely and efficiently.



Cassette holder options

Holders can be produced with single or double sided handles and any depth or shaped rebates to suit all of the unique cassette sizes in use today

Specifications

Overall Size: 2135 x 815mm (965mm wide option)

Patient Surface: 1900 x 680mm (830mm wide option)

Height Range: 430 - 880mm

SWL: 300kg

Castor Size: 150mm twin wheel castors

Package Includes: 100mm premia mattress, folding IV, vertical oxy, foot endboard, emergency CPR backrest release, E-Point steering system, drop down push handles, 2nd handset & USB port

Options: Unique InfineX x-ray system option

Product Codes

Contour Portare-X: MPT-C-PEX

Contour Portare-X Package: MPT-C-PEX-P

150mm Wider Top: MPT-CA-WCT

X-Ray Cassette Holder: MPT-CA-XCH-D

Contour Multi-X

A trauma x-ray stretcher with ability to custom design and build to adapt to any emergency system.



Key Benefits

Smooth and steady: An advanced E-Point steering system and solid, perfectly centred steer wheel ensure the stretcher accurately tracks and turns, while maintaining a strong grip with the floor.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for ED / Imaging

Why Buy?

This stretcher allows our unique InfineX* tracking which gives totally unobstructed head to toe clearance for effortless and efficient x-ray positioning, with an electric backrest that can raise separately to the cassette support surface. We can help you get the most out of your expensive x-ray systems by building to suit the brand and equipment in use. Tailor make the Contour Multi-X with height range choices and shorter, longer, narrower or wider top sizes. Put safety first and select from a range of dropside styles to suit the department requirements.

*Patent pending.

Contour Multi-X



Electric operation

Battery powered touch button control.
Sit and position your patient comfortably and then raise them to the desired work height with ease.



Hydraulic option

Non-powered, manually operated stretcher option.
In remote facilities, or large departments, hydraulic stretchers may be more beneficial, with simple servicing and without the requirement of regular charging or power cords. Includes a tracking digital cassette carriage and holder.



Electric standard

Standard powered stretcher in base configuration.
Sit and position your patient comfortably and then raise them to the desired work height with ease.



Ultra-low height

Easier for patients to get on and off.
The 545mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Excellent high height

Preserve your back and raise a patient up to your best working height.
Lifting the platform height up to 915mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Head and foot tilt

A generous 15 degrees top tilt both directions.
Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Easy customisation

Increase top width, change colours, mattress and dropside styles.
Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.

Specifications

Overall Size: 2135 x 815mm (705mm narrow & 965mm wide options)

Patient Surface: 1900 x 680mm (570mm narrow & 830mm wide options)

Electric Height Range: 545 - 915mm

Hydraulic Height Range: 605 - 930mm

SWL: 300kg

Castor Size: 200mm castors and Ezi-Trak 5th wheel steering

Electric Package Includes: 100mm premia mattress, folding IV, dropdown footboard, vertical oxy, drop down push handles, 2nd handset & USB port

Hydraulic Package Includes: 100mm premia mattress, folding IV, dropdown footboard, vertical oxy & drop down push handles

Product Code

Electric: MPT-C-XE-EB

Electric Package: MPT-C-XE-EB-P

Hydraulic: MPT-C-XH

Hydraulic Package: MPT-C-XH-P

X-Ray Cassette Holder: MPT-CA-XCH-D

Contour Barituff-X

A 500kg SWL extra-wide electric trauma stretcher designed for Bariatric patients.



Key Benefits

Purpose built: Increased componentry strength in the core frame design operates with higher capacity columns to safely lift 500kg without overload.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for ED / Imaging and Bariatric Care

Why Buy?

Having equipment that is robust enough, that you can completely rely on, is the greatest challenge when caring for Bariatric patients. The Contour Barituff-X has increased componentry strength in the core frame design and operates with higher capacity columns to safely lift 500kg without overload.

The stronger and wider top comfortably caters for large patients, with a choice of mattress styles available to suit the type of departmental use. Stretcher stability is assured through an increased base size. Powered drive and weigh scales genuinely reduce manual handling risks to the carer.

This stretcher has been built for Bariatric patients.

Contour Barituff-X



Ultra-low height

Easier for patients to get on and off.
The 600mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Head and foot tilt

A generous 15 degrees top tilt both directions.
Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 980mm (plus an extra 110mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Cassette holder options

Holders can be produced with single or double sided handles and any depth or shaped rebates to suit all of the unique cassette sizes in use today

Specifications

Overall Size: 2135 x 1100mm

Patient Surface: 1900 x 980mm

Mattress: 110mm thick with 150mm thick option

Height Range: 600 - 980mm

SWL: 500kg

Castor Size: 150mm twin wheel castors

Package Includes: 110mm Premia Bariatric mattress, folding IV, weigh scales, vertical oxy, orbit drive, 2nd handset & USB port

Product Codes

Contour Barituff-X: MPT-C-BEX

Contour Barituff-X Package: MPT-C-BEX-P

X-Ray Cassette Holder: MPT-CA-XCH-D

Contour Endosurge

An endoscopy imaging stretcher with a sliding top for quick and easy patient positioning.



Key Benefits

Smooth and steady: An advanced E-Point steering system and solid, perfectly centred steer wheel ensure the stretcher accurately tracks and turns, while maintaining a strong grip with the floor.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for Endoscopy / Colonoscopy & Pain Management

Why Buy?

Take away the effort of re-positioning a heavy C-Arm X-Ray camera. For quick and easy movement of the patient, unlock the stretcher top and slide it a generous distance back or forwards.

Our aim is to provide more freedom and take strain off your backs. The cutaway base cover allows more clearance for the camera under the stretcher whilst bringing it closer to the floor for a lower working height.

Contour Endosurge



Slide to head

Unlock and move the top 140mm towards the head end from centre position.

Release the top lock and easily slide the stretcher top back or forwards to position your patient without having to significantly move the heavy C-arm imaging camera.



Hydraulic option

Non-powered, manually operated stretcher option.

In some cases, such as small or remote facilities, or large high pressure departments, hydraulically operated stretchers may be more beneficial, with simple mechanical servicing and without the requirement of regular charging or power cords.



Slide to foot

Unlock and move the top 260mm towards the foot end from the centre position.

Release the top lock and easily slide the stretcher top back or forwards to position your patient without having to significantly move the heavy C-arm imaging camera.

Specifications

Overall Size: 2135 x 815mm (705mm narrow & 965mm wide options)

Patient Surface: 1900 x 680mm (570mm narrow & 830mm wide options)

Electric Height Ranges: 570 - 940mm

Hydraulic Height Ranges: 630 - 960mm

SWL: 300kg

Castor Size: 200mm castors and Ezi-Trak 5th wheel steering

Sliding Top: 400mm slide; 140mm to head and 260mm to foot from centre

Electric Package Includes: 100mm premia mattress, folding IV, vertical oxy, dropdown footboard, drop down push handles, emergency CPR backrest release, split dropsides, 2nd handset & USB port

Hydraulic Package Includes: 100mm premia mattress, folding IV, vertical oxy, dropdown footboard, drop down push handles & split dropsides

Options: Split dropside option (shown) and mid-position catch upgrades

Product Code

Electric: MPT-C-EE

Electric Package: MPT-C-EE-P

Hydraulic: MPT-C-EH

Hydraulic Package: MPT-C-EH-P

Contour Endo-X

The endoscopy imaging stretcher with a 'C' base for complete freedom of camera movement.



Key Benefits

Optimal design: Increased column spacing with left or right hand chassis

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for Endoscopy & Imaging

Why Buy?

An increased distance between the lifting columns of this stretcher allows excellent coverage with the C-Arm X-Ray camera, avoiding the risks of re-positioning and moving patients during procedures. Complete freedom of movement is also facilitated by the single rear bar joining the stretcher base, which creates ample space for the camera and leaves no restriction to the working height required by staff.

Contour Endo-X



Base option



Ultra-low height

Very safe and easy for patients to get on and off.

The 570mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 940mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Easy customisation

Increase top width, change colours, mattress and dropside styles.

Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.

Specifications

Overall Size: 2135 x 815mm (705mm narrow & 965mm wide options)

Patient Surface: 1900 x 680mm (570mm narrow & 830mm wide options)

Height Range: 570 - 940mm

SWL: 300kg

Castor Size: 200mm large bearing castors

Package Includes: 100mm premia mattress, folding IV, dropdown footboard, vertical oxy, split drop sides, emergency CPR backrest release, C-shaped base for maximum camera clearance, 2nd handset & USB port

Options: Split dropside option (shown) and mid-position catch upgrades

Product Codes

Contour Endo-X: MPT-C-EX

Contour Endo-X Package: MPT-C-EX-P

Split fold dropside option: MPT-CA-SDS

Notes
