

Contour Recline Barituff

The Barituff is a 500kg SWL bariatric chair bed with an extra-wide top and deeper seat.



Key Benefits

Low & safe: The ultra-low height is combined with dropsides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patients legs or clothes.

Purpose built: Increased componentry strength in the core frame design operates with higher capacity columns to safely lift 500kg without overload.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for Day Surgery, ED, Podiatry & General Ward

Why Buy?

Having equipment that is robust enough, that you can completely rely on, is the greatest challenge when caring for bariatric patients.

We have increased componentry strength in the core frame design which operates with higher capacity columns to safely lift 500kg without overload. A wider top with a deeper seat has been standardized to comfortably cater for large patients. Heavy-duty, dark coloured vinyl is used on the mattress cushions to combat wear and staining. This stretcher chair has been built for bariatric patients.

v-20122023

Contour Recline Barituff



Low and safe

Excellent 480mm low height with fold-under drop sides.

The ultra-low height is combined with drop sides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patients legs or clothes.



Easy egress

Forward tilt facilitates getting on and off.

For aged or disabled patients, a gentle incline and the long drop side armrests may assist in safer transfer.



Lies flat

The backrest and leg rest move independently or together from chair to bed.

You now have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Smart transfer

Conveniently stored handles unfold for staff friendly transportation.

Lock in the double steer system, unfold the push handles and move forward effortlessly on the large castors. The unique handle design correctly positions the carer while preventing unwanted forces being applied directly to the patient.



Large top design

The Barituff has a 200mm wider x 140mm deeper seat than the standard Recline.

The seat size is increased to 850mm wide x 600mm deep, easily accommodating larger patients. Use our range of colour coding choices if required to assist in quick identification.



Clear handset

Large simple functions with colour assistance.

Position the patient quickly and easily with the touch of a button - either by independently moving sections or applying one touch contouring.

Specifications

Size: 980mm wide x 1270mm upright (1750mm lying flat) overall

Seat width: 850mm wide x 600mm deep seat

Mattress size: 110mm pressure care mattress cushion with heavy duty cover

Height Range: 480 - 880mm with lie-flat and forward/reverse tilt

SWL: 500kg

Castor size: Twin wheel 150mm castors

Package Includes: Standard IV, vertical oxy holder, drop down footrest & meal tray

Product Code

Contour Recline Barituff: MPC-C-REB

Contour Recline Barituff Package: MPC-C-REB-P

v-20122023