

# Contour Recline Endo

Enhance patient support with ease of access for physicians during procedures such as endoscopy.



## Key Benefits

**Safety & support:** The unique split dropside system allows staff the ability to drop one section without losing critical support for your patient and offers extra length.

**Space saver:** Fit more chairs in the same space and reduce patient transfers, time and OH&S risks.

**Hygiene conscious:** Smooth, gap free platform aids cleaning and helps prevent fluid ingress to meet the highest infection control standards.

**Supercharged:** Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

**Longest warranty:** Peace of mind with a lifetime of local support, easy maintenance and fast guaranteed parts.

## Application

Ideal for Endoscopy / Colonoscopy

## Why Buy?

Make better use of valuable space and streamline your department with the Modsel procedure chairs. They are easy to move in and out of tight areas meaning a higher patient throughput and reduced waiting times. Patients can be positioned easily from sitting upright to a full recline, or lay completely flat with the stability of twin columns under the chair which can be used for CPR or specialised procedures. The unique twin dropside system\* ensures extra safety and efficiency for both the patient and carer. \*International design patents pending

v-02112023

# Contour Recline Endo



## Low and safe

Excellent 530mm low height with fold-under dropsides.

The ultra-low height is combined with dropsides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patients legs or clothes.



## Safer access

Lower the dropsides independently to ensure clear access to the patient without reducing side support.



## Clear handset

Large simple functions with colour assistance.

Position the patient quickly and easily with the touch of a button - either by independently moving sections or applying one touch contouring.



## Lies flat and tilts

The backrest and leg rest move independently or together from chair to bed with trendelenburg.

You now have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat.



## Smart transfer

Conveniently stored handles unfold for staff friendly transportation.

Lock in the double steer system, unfold the push handles and move forward effortlessly on the large castors. The unique handle design correctly positions the carer while preventing unwanted forces being applied directly to the patient.



## Oxy holder options

Choice of backrest or dropside mount oxy bottle mounting.

## Specifications

**Size:** 850 x 1200mm overall (1850mm lying flat)

**Seat width:** Generous patient surface width of 750mm with 100mm thick premium, removable mattress

**Height Range:** 530 - 910mm with lie-flat and forward / reverse tilt

**SWL:** 300kg

**Castor Size:** Twin wheel 150mm castors with double steer

**Package Includes:** Standard IV, vertical oxy holder, drop down footrest, glide-away dropsides with no protrusions & meal tray

## Product Code

Contour Recline Endo: MPC-C-REN

Contour Recline Endo Package: MPC-C-REN-P

v-02112023