

# Contour Recline

A versatile chair to streamline medical examination, procedures and recovery across many healthcare facilities.



## Key Benefits

**Save Space:** Fit more chairs in the same space by increasing throughput with faster recoveries from patients sitting up right.

**Hygiene conscious:** Smooth, gap free platform aids cleaning and helps prevent fluid ingress to meet the highest infection control standards.

**Supercharged:** Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

**Custom quality:** Manufactured in Australia with quality, care and customisation capability.

**Longest warranty:** Peace of mind with a lifetime of local support, easy maintenance and fast guaranteed parts.

## Application

Ideal for Day Surgery, Recovery, Transit Lounge, Fast Track ED & Cardiac

## Why Buy?

Make better use of your valuable space and streamline your department with Modsel's compact procedure chairs. They are easy to move in and out of tight areas meaning a higher patient throughput and reduced waiting times. Patients can be positioned easily from sitting upright to a full recline, or lay completely flat with the stability of twin columns under the chair which can be used for CPR or specialised procedures. We've taken care with every aspect of design to ensure comfort and efficiency for both the patient and carer. Pop-up push handles ensure best operator positioning while preventing transfer movement to the patient.

v-31012023

# Contour Recline



## Low and safe

Excellent 460mm low height with fold-under drop-sides.

The ultra-low height is combined with drop-sides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patients legs or clothes.



## Easy egress

Forward tilt facilitates getting on and off.

For aged or disabled patients, a gentle incline and the long drop-side armrests may assist in safer transfer.



## Lies flat

The backrest and leg rest move independently or together from chair to bed.

You now have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat.



## Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



## Smart transfer

Conveniently stored handles unfold for staff friendly transportation.

Lock in the double steer system, unfold the push handles and move forward effortlessly on the large castors. The unique handle design correctly positions the carer while preventing unwanted forces being applied directly to the patient.



## Wide option

100mm wider top option available.

Change the seat width from 650mm to 750mm wide to easily accommodate larger patients. Use our range of colour coding choices if required to assist in quick identification.



## Clear handset

Large simple functions with colour assistance.

Position the patient quickly and easily with the touch of a button - either by independently moving sections or applying one touch contouring.

## Specifications

**Size:** 780 x 1150mm overall (1750mm lying flat)

**Seat width:** 650mm with 100mm wider option available with 100mm thick premium, removable mattress

**Height Range:** 460 - 840mm with lie-flat and forward / reverse tilt

**SWL:** 300kg

**Castor Size:** Twin wheel 150mm castors with double steer

**Package Includes:** Standard IV, vertical oxy holder, dropdown footrest & meal tray

## Product Codes

Contour Recline: MPC-C-RET

Contour Recline Package MPC-C-RET-P

750mm wide top option: MPC-CA-WCT

v-31012023