

I-move

A height adjustable rounds trolley with storage that is light and easy to move in tight spaces.



Key Benefits

Compact & simple: Light, easy to move and unobtrusive in tight spaces is the perfect partner to make sure your nursing rounds are efficient and easy.

Tech support: Upgrade your I-move to support any and all of your IT requirements required for the job.

Safe & secure: Option to add a Kensington lock to restrain valuable laptop hardware to keep safe as well as cable management options.

Why Buy?

The i-move is an affordable and height adjustable rounds trolley. When it comes to moving laptops, tablets or simply pens and paper along with other clinical items, the only practical solution has been the traditional metal dressing trolley.

How good would it be to have a trolley that is light and easy to move, unobtrusive in tight spaces and that you can comfortably sit or stand at? Coupled with this is a convenient storage shelf and our ability to customise the top size and shape to totally suit your requirements!

v-07052024

I-move



Storage space

Convenient storage space with retaining lip across the front.



Customisable

Use our standard work surface or customise to suit your specific requirements.



Height adjustable

Simply squeeze the lever under the top and effortlessly adjust to your desired work height.



Storage compartment

Add a removable stainless steel storage compartment for cables / pens etc.



Power supply holder

Securely hold a laptop power supply for easy charging.



Kensington lock

Restrain valuable laptop hardware to keep it safe.



Add your IT

Add your own holders to suit the specific IT hardware in use.

Specifications

Size: Compact, easy to clean top - 550mm wide x 560mm deep

Height adjustable: From 840mm to 1140mm

Storage size: Convenient 100mm storage space with lip

Product Codes

I-Move Trolley; MCT-IM-LHA

Removable Storage Container - MCT-IMA-SC

Power Supply Holder - MCT-IMA-PCH

Kensington Lock - MCT-IMA-KL

v-07052024