

# Contour Barituff

A 500kg SWL extra-wide electric stretcher designed for Bariatric patients.



## Key Benefits

**Purpose built:** Increased componentry strength in the core frame design operates with higher capacity columns to safely lift 500kg without overload.

**Effortless control:** Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

**Supercharged:** Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

**Custom quality:** Manufactured in Australia with speed, quality, and care, plus customisation capability.

**Longest warranty:** Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

## Application

Ideal for Day Surgery, ED, General Ward & Bariatric Care

## Why Buy?

Having equipment that is robust enough, that you can completely rely on, is the greatest challenge when caring for Bariatric patients. The Contour Barituff has increased componentry strength in the core frame design and operates with higher capacity columns to safely lift 500kg without overload.

The stronger and wider top comfortably caters for large patients, with a choice of mattress styles available to suit the type of departmental use. Stretcher stability is assured through an increased base size. Powered drive and weigh scales genuinely reduce manual handling risks to the carer. This stretcher has been built for Bariatric patients.

v-11122023

# Contour Barituff



## Ultra-low height

Easier for patients to get on and off.  
The 470mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



## Head and foot tilt

A generous 15 degrees top tilt both directions.  
Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



## Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 850mm (plus an extra 110mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



## Easy customisation

Increase top width, change colours, mattress and dropside styles.

Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.

## Specifications

**Overall Size:** 2135 x 1100mm

**Patient Surface:** 1900 x 980mm

**Mattress:** 110mm thick with 150mm thick option

**Height Range:** 470 - 850mm

**SWL:** 500kg

**Castor Size:** 150mm twin wheel castors

**Package Includes:** 110mm premia bari mattress, folding IV, vertical oxy, weigh scales, orbit drive, emergency CPR backrest release, 2nd handset & USB port

## Product Code

Contour Barituff: MPT-C-BE

Contour Barituff Package: MPT-C-BE-P

v-11122023