

Contour Multi-X

A trauma x-ray stretcher with ability to custom design and build to adapt to any emergency system.



Key Benefits

Smooth and steady: An advanced E-Point steering system and solid, perfectly centred steer wheel ensure the stretcher accurately tracks and turns, while maintaining a strong grip with the floor.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for ED / Imaging

Why Buy?

This stretcher allows our unique InfineX* tracking which gives totally unobstructed head to toe clearance for effortless and efficient x-ray positioning, with an electric backrest that can raise separately to the cassette support surface. We can help you get the most out of your expensive x-ray systems by building to suit the brand and equipment in use. Tailor make the Contour Multi-X with height range choices and shorter, longer, narrower or wider top sizes. Put safety first and select from a range of dropside styles to suit the department requirements.

*Patent pending.

v-08122023

Contour Multi-X



Electric operation

Battery powered touch button control.
Sit and position your patient comfortably and then raise them to the desired work height with ease.



Hydraulic option

Non-powered, manually operated stretcher option.
In remote facilities, or large departments, hydraulic stretchers may be more beneficial, with simple servicing and without the requirement of regular charging or power cords. Includes a tracking digital cassette carriage and holder.



Electric standard

Standard powered stretcher in base configuration.
Sit and position your patient comfortably and then raise them to the desired work height with ease.



Ultra-low height

Easier for patients to get on and off.
The 545mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Excellent high height

Preserve your back and raise a patient up to your best working height.
Lifting the platform height up to 915mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Head and foot tilt

A generous 15 degrees top tilt both directions.
Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Easy customisation

Increase top width, change colours, mattress and dropside styles.
Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.

v-08122023

Specifications

Overall Size: 2135 x 815mm (705mm narrow & 965mm wide options)

Patient Surface: 1900 x 680mm (570mm narrow & 830mm wide options)

Electric Height Range: 545 - 915mm

Hydraulic Height Range: 605 - 930mm

SWL: 300kg

Castor Size: 200mm castors and Ezi-Trak 5th wheel steering

Electric Package Includes: 100mm premia mattress, folding IV, dropdown footboard, vertical oxy, drop down push handles, 2nd handset & USB port

Hydraulic Package Includes: 100mm premia mattress, folding IV, dropdown footboard, vertical oxy & drop down push handles

Product Code

Electric: MPT-C-XE-EB

Electric Package: MPT-C-XE-EB-P

Hydraulic: MPT-C-XH

Hydraulic Package: MPT-C-XH-P

X-Ray Cassette Holder: MPT-CA-XCH-D

v-08122023