

Contour Portare-X

An ultra-low electric trauma stretcher with totally unobstructed x-ray clearance.



Key Benefits

Unobstructed x-rays: Supported by just four outside corners, our unique x-ray system gives total freedom of movement for cassettes - even with an electric backrest!

Safe & low: The 430mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting.

Smooth & steady: Operated by a very light touch pedal action, this wheel positively pre-loads onto the floor to ensure strong grip and rides up and down effortlessly over bumps and hollows to keep traction.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for ED & Imaging

Why Buy?

Helping patients on and off a trauma stretcher is safer and easier when the bed can be lowered down to 430mm, around the standard office chair height! This stretcher can feature our unique InfineX* system which allows totally unobstructed head to toe clearance for effortless and efficient x-ray cassette positioning. Smoothly moving patients helps improve their experience and our E-Point steering delivers a strong grip to the floor with a solid, perfectly centered wheel to allow the stretcher to always accurately track and turn. The Modsel x-ray stretchers are designed to be a reliable, cost-effective performer, no matter what the shift brings.

*Patent pending.

v-21122023

Contour Portare-X



Unobstructed x-ray tunnel

100% x-ray clearance of the entire body. Supported by just four outside corners, our unique patented InfineX system gives total freedom of movement for x-ray cassettes - even with an electric backrest! There are two styles available, another that allows the backrest to raise separately to the lower cassette support surface.



Excellent high height

Quickly raise a patient up to the best working height. Lifting the platform height up to 880mm helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



E-Point steering

A perfectly centred steer wheel ensures the stretcher always accurately tracks and turns. Operated by a very light touch pedal action, this wheel positively pre-loads onto the floor to ensure strong grip and rides up and down effortlessly over bumps and hollows to keep traction. Turning left or right from either end gives the same steering result.



Cassette holder options

Holders can be produced with single or double sided handles and any depth or shaped rebates to suit all of the unique cassette sizes in use today



Ultra-low height

Safer and easier for unwell patients to get on and off. The 430mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Head and foot tilt

A generous 15 degrees top tilt both directions. Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Easy customisation

Increase top width, change colours, mattress and dropside styles. Every trauma department has their own unique requirements to be able to function safely and efficiently.

v-21122023

Specifications

Overall Size: 2135 x 815mm (965mm wide option)

Patient Surface: 1900 x 680mm (830mm wide option)

Height Range: 430 - 880mm

SWL: 300kg

Castor Size: 150mm twin wheel castors

Package Includes: 100mm premia mattress, folding IV, vertical oxy, foot endboard, emergency CPR backrest release, E-Point steering system, drop down push handles, 2nd handset & USB port

Options: Unique InfineX x-ray system option

Product Codes

Contour Portare-X: MPT-C-PEX

Contour Portare-X Package: MPT-C-PEX-P

150mm Wider Top: MPT-CA-WCT

X-Ray Cassette Holder: MPT-CA-XCH-D

v-21122023