

Contour Portare

An all-purpose electric transport and recovery stretcher with excellent low height.



Key Benefits

Safer solution: The 400mm low platform height of this stretcher takes away the requirement to use step stools or unnecessary heavy lifting.

Effortless control: Sit and position your patient comfortably and then raise and lower them to the desired work height with the touch of a button.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Transfer stretcher suitable for most medical departments

Why Buy?

Helping patients on and off a general purpose stretcher has never been this easy or safe with a design that allows the mattress base down to 400mm, lower than a standard office chair height! Transport must always be considered and with our new E-Point system, a strong grip with the floor and solid, perfectly centered steer wheel ensures the stretcher always accurately tracks and turns. We bring Modsel stretcher design into department areas such as emergency, theatre and day surgery helping to maximise the number of working units available.

v-01032023

Contour Portare



Ultra-low height

Very safe and easy for patients to get on and off.

The 400mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



E-Point steering

A perfectly centred steer wheel ensures the stretcher always accurately tracks and turns.

Operated by a very light touch pedal action, this wheel positively pre-loads onto the floor to ensure strong grip and rides up and down effortlessly over bumps and hollows to keep traction. Turning left or right from either end gives the same steering result.



Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 830mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Easy customisation

Increase top width, change colours, mattress and dropside styles.

Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.

Specifications

Overall Size: 2135 x 815mm (965mm wide option)

Patient Surface: 1900 x 680mm (830mm wide option)

Height Range: 400 - 830mm

SWL: 300kg

Castor Size: 150mm twin wheel castors

Package Includes: 100mm premia mattress, folding IV, vertical oxy, foot endboard, drop down push handles, emergency CPR backrest release, 2nd handset & USB port

Product Code

Contour Portare: MPT-C-PE2

Contour Portare Package: MPT-C-PE2-P

150mm Wider Top: MPT-CA-WCT

v-01032023