

# Contour Jazz

An Australian made transport and recovery stretcher that's a joy to use and designed for supreme patient safety and comfort.



## Key Benefits

**Smooth and steady:** An advanced E-Point steering system and solid, perfectly centred steer wheel ensure the stretcher accurately tracks and turns, while maintaining a strong grip with the floor

**Safety and support:** The ultra-low 400mm height makes it safer and easier for patients to get on and off – no need for step stools or unnecessary heavy lifting.

**Effortless control:** Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

**Bright and light:** Available in a choice of four fun and colourful themes to brighten their day (and yours).

**Longest warranty:** Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

## Application

Ideal for: Emergency, Theatre & Day Surgery

## Why Buy?

Keeping patients safe, secure, and comfortable has never been easier with the Contour Jazz transport and recovery stretcher. With full, touch button electronic control and innovative E-Point steering system, the stretcher offers a smooth ride for patients and the kind of ease of operation and manoeuvrability that every caring professional craves. The Contour Jazz has been expertly designed to ensure you'll never miss a beat in even the most challenging environment or circumstances. Featuring a fast-action CPR backrest release, strong grip with the floor and an ultra-low height range that allows the mattress base down to 400mm. That's lower than a standard office chair height! Available in a range of colours and themes.

v-20122023

# Contour Jazz



Ready for a range of accessories



Ultra-low height

Very safe and easy for patients to get on and off.

The 400mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Music Theme



Sea Theme



Jungle Theme



Space Theme



Fixed Control Panel

A fixed control panel on the foot end of the stretcher prevents unwanted misuse.

v-20122023

## Specifications

**Overall Size:** 2135 x 815mm (965mm wide option)

**Patient Surface:** 1900 x 680mm (830mm wide option)

**Height Range:** 400mm-830mm

**SWL:** 300kg

**Castor Size:** 150mm twin wheel castors

**Package Includes:** Handset and end controls, 100mm Premia mattress, emergency CPR back rest release, foot end board, drop down push handles, vertical oxy holder, folding IV, 16 riser dropsides and USB Port.

## Product Code

MPT-C-JAZZ

## Decal Codes

MPT-CJA-MUSIC

MPT-CJA-SEA

MPT-CJA-JUNGLE

MPT-CJA-SPACE

v-20122023