



General Purpose Stretchers



General Purpose Stretchers

We are medical bed suppliers building patient transport stretchers that not only cater to your department's operational efficiency, but also to your workload and safety. Our stretchers are designed to last and are a long-term, safe investment.

The trusted option for medical beds

Our complete range of stretchers and beds are all designed to meet the demanding needs of healthcare environments. Get in touch with Modsel on 1300 785 790 to have any questions answered.



Contour Recline



Contour Energy



Contour Portare



Contour Conceal



Contour Jazz




Contour Classic



Contour Exami



Contour Deluxe



Stretcher Mattress



Orbit-Drive



Oxy Holders



Dropdown Push Handles



End Rail



Dropdown Foot Board



Change Table

Contour Recline

A versatile chair to streamline medical examination, procedures and recovery across many healthcare facilities.



Key Benefits

Save Space: Fit more chairs in the same space by increasing throughput with faster recoveries from patients sitting up right.

Hygiene conscious: Smooth, gap free platform aids cleaning and helps prevent fluid ingress to meet the highest infection control standards.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with quality, care and customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance and fast guaranteed parts.

Application

Ideal for Day Surgery, Recovery, Transit Lounge, Fast Track ED & Cardiac

Why Buy?

Make better use of your valuable space and streamline your department with Modsel's compact procedure chairs. They are easy to move in and out of tight areas meaning a higher patient throughput and reduced waiting times. Patients can be positioned easily from sitting upright to a full recline, or lay completely flat with the stability of twin columns under the chair which can be used for CPR or specialised procedures. We've taken care with every aspect of design to ensure comfort and efficiency for both the patient and carer. Pop-up push handles ensure best operator positioning while preventing transfer movement to the patient.

Contour Recline



Low and safe

Excellent 460mm low height with fold-under drop-sides.

The ultra-low height is combined with drop-sides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patient's legs or clothes.



Easy egress

Forward tilt facilitates getting on and off.

For aged or disabled patients, a gentle incline and the long drop-side armrests may assist in safer transfer.



Lies flat

The backrest and leg rest move independently or together from chair to bed.

You now have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Smart transfer

Conveniently stored handles unfold for staff friendly transportation.

Lock in the double steer system, unfold the push handles and move forward effortlessly on the large castors. The unique handle design correctly positions the carer while preventing unwanted forces being applied directly to the patient.



Wide option

100mm wider top option available.

Change the seat width from 650mm to 750mm wide to easily accommodate larger patients. Use our range of colour coding choices if required to assist in quick identification.



Clear handset

Large simple functions with colour assistance.

Position the patient quickly and easily with the touch of a button - either by independently moving sections or applying one touch contouring.

Specifications

Size: 780 x 1150mm overall (1750mm lying flat)

Seat width: 650mm with 100mm wider option available with 100mm thick premium, removable mattress

Height Range: 460 - 840mm with lie-flat and forward / reverse tilt

SWL: 300kg

Castor Size: Twin wheel 150mm castors with double steer

Package Includes: Standard IV, vertical oxy holder, dropdown footrest & meal tray

Product Codes

Contour Recline: MPC-C-RET

Contour Recline Package MPC-C-RET-P

750mm wide top option: MPC-CA-WCT

Contour Energy

A transport chair with ZERO effort, 4-way power drive reducing the risk of strain during patient transport.



Key Benefits

Ergonomics & intuitive: Drive effortlessly in any direction with only the slightest of touch, paired with perfect positioning of the drive controls for ease of use.

Easy access: The large footplates offer generous support to your patients and lift quickly for clear egress.

Powerful movement: A rotating electric drive wheel provides excellent traction, while center steering makes patient transport easier over longer distances.

Why Buy?

Take all the human effort out of pushing and maneuvering chairs with an intuitive and intelligent drive system designed for easy use by any operator. Natural movement controls the speed and direction of the chair with the ability to stop very quickly by pulling back. Help prevent the most common carer strain injuries by taking the back twist out of moving heavy weight forwards, backwards and sideways. Patients are positioned comfortably with a gentle recline and feet fully supported. We've taken care with every aspect of design to ensure comfort and efficiency for both the patient and carer. Customise for a wide range of healthcare areas with relevant accessories.

Contour Energy



Ergonomic

Drive effortlessly in any direction by only light touch.



Intuitive

Perfect positioning of the drive controls for the operator.



Support

The large foot plates offer generous support to the patient.



Easy access

The large foot plates lift quickly and easily for clear egress.



Full clearance

Complete frontal clearance by lifting and rotating foot plates outwards.



Removable armrests

The side armrests can be fully removed for multi direction egress.



Versatile

Use the storage shelf during transport or customise with accessories such as IV poles, oxy holders and seat belt.

Specifications

Size: 750mm W x 1300mm L x 1180mm H

Seat Size: 650mm W x 520mm D x 550mm H

SWL: 300kg

Castor Size: 150mm

Package Includes: standard IV & vertical oxy holder

Product Code

Contour Energy: MPC-C-ENT-EDU

Contour Energy Package: MPC-C-ENT-EDU-P

Contour Portare

An all-purpose electric transport and recovery stretcher with excellent low height.



Key Benefits

Safer solution: The 400mm low platform height of this stretcher takes away the requirement to use step stools or unnecessary heavy lifting.

Effortless control: Sit and position your patient comfortably and then raise and lower them to the desired work height with the touch of a button.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Transfer stretcher suitable for most medical departments

Why Buy?

Helping patients on and off a general purpose stretcher has never been this easy or safe with a design that allows the mattress base down to 400mm, lower than a standard office chair height! Transport must always be considered and with our new E-Point system, a strong grip with the floor and solid, perfectly centered steer wheel ensures the stretcher always accurately tracks and turns. We bring Modsel stretcher design into department areas such as emergency, theatre and day surgery helping to maximise the number of working units available.

Contour Portare



Ultra-low height

Very safe and easy for patients to get on and off.

The 400mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



E-Point steering

A perfectly centred steer wheel ensures the stretcher always accurately tracks and turns.

Operated by a very light touch pedal action, this wheel positively pre-loads onto the floor to ensure strong grip and rides up and down effortlessly over bumps and hollows to keep traction. Turning left or right from either end gives the same steering result.



Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 830mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Easy customisation

Increase top width, change colours, mattress and dropside styles.

Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.

Specifications

Overall Size: 2135 x 815mm (965mm wide option)

Patient Surface: 1900 x 680mm (830mm wide option)

Height Range: 400 - 830mm

SWL: 300kg

Castor Size: 150mm twin wheel castors

Package Includes: 100mm premia mattress, folding IV, vertical oxy, foot endboard, drop down push handles, emergency CPR backrest release, 2nd handset & USB port

Product Code

Contour Portare: MPT-C-PE2

Contour Portare Package: MPT-C-PE2-P

150mm Wider Top: MPT-CA-WCT

Contour Conceal

A purpose-built concealment stretcher for mortuary transport with superior functionality and customisable features.



Key Benefits

Careful concealment: Fully customisable concealment stretcher with high quality, perfectly fitted vinyl cover in subtle grey or printed in one of three restful, full colour art designs.

Effortless control: Easy-drop flushfold sides and adjustable stretcher height ensures ease of transfer onto the Contour Conceal. Electric and hydraulic options available.

Smooth and steady: Robust, two-column design with 200mm castors and Ezi-Trak 5th wheel steering ensures stretcher stability and manoeuvrability in-transit

Easy clean: The Contour Conceal's smooth, one-piece platform and removeable, washable vinyl cover combine to make cleaning and hygiene control simple and effective.

Application

Ideal for Wards, Emergency, Theatre, Intensive Care & Palliative

Why Buy?

The Contour Conceal has been sensitively designed to assist you to provide patients who have passed away with a safe, secure, and dignified passage from wards to mortuaries. Offering full concealment, with your choice of subtle grey or calming rainforest, mountain or lakeside theme covers, the platform stretcher ensures each transfer is carried out with the utmost of ease, efficiency and respect. The twin column design is available as either an electric or hydraulic stretcher with a narrow or wide platform top. Both options are able to be adjusted to several height ranges and feature flushfold dropsides to facilitate the patient's smooth transfer onto the stretcher.

Contour Conceal



Rainforest Theme



Lake Theme



Mountain Theme



Clear access

Lift off the heavy duty washable vinyl cover with full length access from both sides. Store the cover supports on convenient hooks.



Custom Options

Upgrade your Contour Conceal with orbit-drive, out door base or your choice of a grey or printed theme cover



Top width options

Narrow and wider top options available



Hydraulic option

Choose manual height adjustment if electric operation is not required

Specifications

Overall Size: 2135 x 815mm (705mm narrow & 965mm wide options)

Patient Surface: 1900 x 700mm (590mm narrow & 850mm wide options)

Electric Height Range: 440-810mm (standard) up to 540-910mm

Hydraulic Height Range: 525-850mm

SWL: 300kg

Castor Size: 200mm castors & Ezi-Trak 5th wheel steering

Options: Power drive system and outdoor base upgrades available

Product Code

Electric - MPT-C-CCE

Hydraulic - MPT-C-CCH

Printed Cover Codes

MPT-CCA-LAKE

MPT-CCA-MOUNT

MPT-CCA-RAIN

Contour Jazz

An Australian made transport and recovery stretcher that's a joy to use and designed for supreme patient safety and comfort.



Key Benefits

Smooth and steady: An advanced E-Point steering system and solid, perfectly centred steer wheel ensure the stretcher accurately tracks and turns, while maintaining a strong grip with the floor

Safety and support: The ultra-low 400mm height makes it safer and easier for patients to get on and off – no need for step stools or unnecessary heavy lifting.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Bright and light: Available in a choice of four fun and colourful themes to brighten their day (and yours).

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for: Emergency, Theatre & Day Surgery

Why Buy?

Keeping patients safe, secure, and comfortable has never been easier with the Contour Jazz transport and recovery stretcher. With full, touch button electronic control and innovative E-Point steering system, the stretcher offers a smooth ride for patients and the kind of ease of operation and manoeuvrability that every caring professional craves. The Contour Jazz has been expertly designed to ensure you'll never miss a beat in even the most challenging environment or circumstances. Featuring a fast-action CPR backrest release, strong grip with the floor and an ultra-low height range that allows the mattress base down to 400mm. That's lower than a standard office chair height! Available in a range of colours and themes.

Contour Jazz



Ready for a range of accessories



Ultra-low height

Very safe and easy for patients to get on and off.

The 400mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Music Theme




Sea Theme



Jungle Theme



Space Theme



Fixed Control Panel

A fixed control panel on the foot end of the stretcher prevents unwanted misuse.

Specifications

Overall Size: 2135 x 815mm (965mm wide option)

Patient Surface: 1900 x 680mm (830mm wide option)

Height Range: 400mm-830mm

SWL: 300kg

Castor Size: 150mm twin wheel castors

Package Includes: Handset and end controls, 100mm Premia mattress, emergency CPR back rest release, foot end board, drop down push handles, vertical oxy holder, folding IV, 16 riser dropsides and USB Port.

Product Code

MPT-C-JAZZ

Decal Codes

MPT-CJA-MUSIC

MPT-CJA-SEA

MPT-CJA-JUNGLE

MPT-CJA-SPACE

Contour Classic

The most flexible custom design platform available for a transport and recovery stretcher.



Key Benefits

Customisable design: Choose from a variety of dropside options and electric or hydraulic versions available.

Effortless control: Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for Day Surgery, Recovery, Transit Lounge & Cardiac

Why Buy?

The most flexible platform of all our stretcher range. The twin column design allows for several height range choices - we can build in a low height range for patient egress care or higher height ranges for better equipment or surgeon access with less strain during procedures.

Tailor make the Contour Classic with shorter, longer, narrower and wider top sizes and select from a range of dropside styles to suit specialised procedures or the patient type you are working with. To simply put it, there is a high level of customisation possible with many accessory mounting options and solutions.

Contour Classic



Hydraulic option

Non-powered, manually operated stretcher option.

In some cases, such as small or remote facilities, or large high pressure departments, hydraulically operated stretchers may be more beneficial, with simple mechanical servicing and without the requirement of regular charging or power cords.



Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 810mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Easy customisation

Increase top width, change colours, mattress and dropside styles.

Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.



Ultra-low height

Very safe and easy for patients to get on and off.

The 440mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.

Specifications

Overall Size: 2135 x 815mm (705mm narrow & 965mm wide options)

Patient Surface: 1900 x 680mm (570mm narrow & 830mm wide options)

Electric Height Ranges: 440 - 810 (std) up to 540 - 910mm

Hydraulic Height Range: 525 - 850mm

SWL: 300kg

Castor Size: 200mm castors and Ezi-Trak 5th wheel steering

Electric Package Includes: 100mm premia mattress, folding IV, vertical oxy, dropdown footboard, drop down push handles, 2nd handset & USB port

Hydraulic Package Includes: 100mm premia mattress, folding IV, vertical oxy, dropdown footboard & drop down push handles

Product Code

Electric: MPT-C-CE

Electric Package: MPT-C-CE-P

Hydraulic: MPT-C-CH

Hydraulic Package: MPT-C-CH-P

Contour Exami

An economical, customisable, easy to move examination table available with or without dropsides.



Key Benefits

Safe & low: The 400mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting.

Excellent height: Electrically adjust the patient height up to 930mm helps to examine or treat your patient without having to stoop over, ensuring easy access and less stress.

Comfortable: 100mm fluid proof Protran mattress with a welded, pressure relieving cover.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for medical clinics, examination rooms, schools, prisons & mining sites

Why Buy?

Helping patients on and off an examination table is easy and safe with a design that allows the mattress platform down to 400mm, lower than a standard office chair height! Uniquely, these stretcher examination tables are designed to be easily moved around and re-positioned due to their light-weight build and large castors. There are a range of accessories and customisations available to ensure the unit is designed to be priced and configured best to practically fit into your modern facility.

Contour Exami



Without dropsides

If you require a patient examination table, then this is worth considering. The standard Contour Exami without dropsides gives all the same functionality as an examination table, however has the additional value of colour choices and accessory mounting, while being very easy to move and reposition as required.



Excellent high height

Don't strain your back, raise a patient up to your best working height. Lifting the platform height up to 830mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release option available.



Customise it your way

Upgrades such as wider top, dropsides, base cover, colours and accessories. This is a very flexible platform with excellent customisation available. Change the standard Contour Exami from a basic examination table into an economy patient stretcher.



Ultra-low height

Make it easier and safer for patients to get on and off. The 400mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Head and foot tilt

A generous 15 degrees top tilt both directions. Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Go anywhere

Change from the standard 125mm twin to 175mm single castors for situations such as outdoor. If there is a requirement to wheel over rougher concrete, sand or grass terrains, this option gives larger diameter rubber wheels with 50mm extra ground clearance under the stretcher.

Specifications

Overall Size: 2135 x 815mm (965mm wide option)

Patient Surface: 1900 x 680mm (830mm wide option)

Height Range: 400 to 830mm

SWL: 300kg

Castor Size: 125mm twin wheel castors, 4 brakes and 1 steer

Includes: 100mm Protran mattress, electric backrest that raises to 85 degrees, standard main power electric operation

Options: dropsides, base cover, battery back-up with cable retractor, larger castors and 6 endrail colours

Product Code

Contour Exami (without dropsides); MPT-C-EXAMI

150mm wider top; MPT-CXA-WT

Dropsides; MPT-CXA-5RD

Base cover with storage; MPT-CXA-BC

Battery backup; MPT-CXA-BB

CPR backrest release; MPT-CXA-BQR

175mm large castor option; MPT-CXA-175

Contour Deluxe

A transport stretcher featuring glide away siderails with zero transfer gap.



Key Benefits

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for Day Surgery, Recovery, Transit Lounge & Cardiac

Why Buy?

Safely transfer patients with a zero transfer gap created by siderails that glide down under the top, leaving no obstruction on the side of the stretcher. Gas assistance ensures they are always light and easy for the operator and there is quick access with release catches on both ends.

Provide a greater level of support with these extra length siderails featuring horizontal bars and reduced gaps. Built for a lifetime, the fully welded and polished stainless steel construction leaves no risk of coatings to chip and rust, creating infection hazards or strength reduction.

Contour Deluxe



Hydraulic option

Non-powered, manually operated stretcher option.

In some cases, such as small or remote facilities, or large high pressure departments, hydraulically operated stretchers may be more beneficial, with simple mechanical servicing and without the requirement of regular charging or power cords.



Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 940mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Easy customisation

Increase top width, change colours, mattress and dropside styles.

Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.



Ultra-low height

Easier for patients to get on and off.

The 570mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.

Specifications

Overall Size: 2135 x 815mm (965mm wide option)

Patient Surface: 1900 x 680mm (830mm wide option)

Electric Height Range: 570 - 940mm

Hydraulic Height Range: 630 - 960mm

SWL: 300kg

Castor Size: 200mm castors and Ezi-Trak 5th wheel steering

Includes: Full length gas assisted siderails that are fully welded stainless steel, emergency CPR backrest release, X-ray subtop option

Product Code

Electric - MPT-C-DE

Hydraulic - MPT-C-DH

Stretcher Mattress

High grade mattress with excellent durability and pressure relieving comfort for patient stretchers.

Key Benefits

- Thicknesses that aid egress and imaging
- Protran mattress - single foam layer
- Premia mattress - 3 effective foam layers
- Premia Barituff mattress - 4 support layers (110mm thick 60 - 350kg, 150mm thick 45 - 400kg therapeutic support loads)
- Welded 'double coated' stretch cover
- Non-slip base design
- Permanent care label
- New technology fluid proof zip (excludes 75mm mattress)
- MADE IN AUSTRALIA

Product Code

- 100mm Protran Mattress; MPT-CM-100-PTP
- 100mm Premia Mattress; MPT-CM-100-PRP
- 110mm Barituff Mattress; MPT-CM-110-PRB
- 150mm Barituff Mattress; MPT-CM-150-PRB

Mattress Size



- Standard; 680 W x 1900mm L
- Narrow; 570mm W x 1900mm L
- Wide; 830mm W x 1900mm L
- Barituff; 980mm W x 1900mm L

Mattress Shape

- Square Ends (electric stretcher)
- Truncated End (hydraulic stretcher)

Why Buy

The Contour Premia mattress has been designed to optimize patient weight capacity and comfort by using a selection of quality foam layers to help in both preventing patients from bottoming and pressure reduction. The Contour Protran mattress is recommended for short stay patients with a single medium density foam section to facilitate special patient transfers or clear imaging.



Pressure care design




Welded cover




Sealed zip



Care label



Pressure mapped



Protran

Orbit-Drive

Taking patient transport to a new level in carer welfare with zero effort, 4-way power drive.



Key Benefits

- Forward, back and sideways powered drive
- Intuitive height adjustable controls with emergency stop
- Powerful electric drive centre steer wheel

Orbit-Drive



Adjustable controls

Lower the drive controls down below the mattress height for enhanced patient access and specialised accessory use.



Adaptable

Install the power drive system on many of the Contour stretcher models



Powerful drive

An electric driven wheel pressured onto the floor that provides excellent traction and center steering.



Sideways movement

The unique rotating design allows the stretcher to be driven in any direction - forwards, backwards and sideways!

Why Buy

Take all the human effort out of pushing and jostling stretchers with an intuitive and intelligent drive system designed for easy use by any operator. Natural movement controls the speed and direction of the stretcher with the ability to stop very quickly by pulling back.

Help prevent the most common carer strain injuries by taking the back twist out of moving a stretcher forwards, backwards and sideways. Importantly, we haven't overlooked safety features such as emergency stop and system lock out, which are standard. Add this system to our traditional vertical column stretcher designs such as Contour Barituff, Classic, Deluxe and Multi-X to improve both standard and Bariatric patient transport.

Product Code

- MPT-CA-EDU

Oxy Holders

A selection of coated stainless steel oxy bottle holders that can be mounted vertically or horizontally.

Key Benefits

- Tough, powder coated stainless steel basket construction
- Vertical basket size for up to 120mm diameter bottle; 320mm deep x 150mm wide
- A range of vertical mount positions - low, high with swivel or 45 and 90 degree fixed angles
- Horizontal basket size for up to 140mm diameter bottle; 600mm long x 180mm wide
- Horizontal basket can vary in hanging height depending on stretcher models; typically 280mm

Product Code

- Vertical Rail Mount - MPT-CA-VO2
- Vertical Corner Mount Swivel - MPT-CA-VCO2
- Vertical Corner Mount Fixed 45 - MPT-CA-VCO2-45
- Vertical Corner Mount Fixed 90 - MPT-CA-VCO2-90
- Vertical Corner High Mount Swivel - MPT-CA-VCUO2
- Vertical Corner High Mount Fixed 45 - MPT-CA-VCUO2-45
- Vertical Corner High Mount Fixed 90 - MPT-CA-VCUO2-90
- Horizontal Oxy Holder - MPT-CA-HO2

Why Buy

Due to so many wide and varied departmental requirements for specific oxy bottle sizes and positioning, we have developed a range of options. Larger bottles can be laid horizontally to ensure they are kept out of the way and smaller bottles are commonly mounted vertically. Some prefer 'fixed', non-moving positions and others require 'swinging' options that allow the bottles to be moved aside if they become an obstruction during transport. There is also the option to support bottles vertically both above and below the patient surface.



Vertical Corner Mount



Vertical Corner High Mount



Horizontal Oxy Holder



Horizontal Oxy Holder

Dropdown Push Handles

Strong fold-away handles with rubber grip for pushing the trolley.

Key Benefits

- Strong stainless steel construction
- Lift up and fold down
- Tough hygienic hand grips
- Tall head end handle option - 400mm high x 600mm wide (350mm between handles)
- Foot end handle option - 160mm high x 600mm wide (350mm between handles)

Product Codes

- Tall head end handles; MPT-CA-TDPH
- Foot end handles; PT-CA-DPH

Why Buy

To allow space for the many removable and common accessories used at the foot end of the stretcher, the handles are shorter to fold completely away with no obstruction to the accessory mounting positions. There is an option for a tall version of these handles to be mounted at the head end of the stretcher which provide high pushing points and ample clearance from the moving backrest.



Head end handles folded



Foot end handles



Foot end handles folded

End Rail

An end rail that can help to retain the mattress, patient or be used as a push handle.

Key Benefits

- Strong, 25mm diameter coated stainless steel rail
- 600mm wide x 200mm high
- Fits into foot end or optional head end accessory sockets

Product Codes

- End rail; MPT-CA-FER
- Head end accessory sockets; MPT-CA-ASH (pair)

Why Buy

Quickly drop this rail into the accessory sockets (optional accessory on the head end of a stretcher, excludes Contour Vertex) and use it as a higher strong point to move the stretcher. This rail has a variety of additional uses such retaining a mattress or pillow from sliding off or a bar for a patients feet to push against. These rails can be used on either end of the stretcher or both ends together.



Head end option



Both ends

Dropdown Foot Board

Padded fold-down foot board to increase the length of the patient surface.

Key Benefits

- High density foam covered in welded stretch material
- Tough board construction in a one-piece hospital grade plastic
- Stainless steel lift and lock folding hinge
- 530mm wide x 330mm high (250mm high cushion)
- Adds an extra 310mm to the length of a mattress

Product Code

- MPT-CA-DDF

Why Buy

Padded for patient comfort, the removable dropdown foot board has molded handles for pushing the stretcher and with a simple 'lift and lay' action, folds down to quickly increase the length of the mattress area.



Fold down

Change Table

A purpose built, non-slip change mat on a sturdy, height adjustable and safe patient platform.



Key Benefits

Safe & secure: The fold up dropsides and removable end rails ensure high levels of patient safety when required.

Customisable design: Tailor to your requirements with options such as electric high/low, 300kg SWL and wide tops.

Why Buy?

This change mat expands the use of the Aquarius trolley into a very safe, height adjustable change table with multi-position, removable end rails and full length dropsides that fold back underneath the trolley top for ample work clearance. The mat securely fastens underneath the trolley top to prevent sliding and movement. This trolley is designed for harsh use and able to be rigorously cleaned.

Change Table



Safe design

The fold up dropsides and removable end rails ensure a high level of patient safety when required.



Customise

Tailor to your requirements with options such as electric high/low, 300kg SWL and wide tops.

Specifications

Size: 795 W x 1950mm including bumpers and 25mm high density foam core

Patient Surface Size: 600 W x 1820mm Long

Height Range: 60 to 900mm (depending on trolley model)

SWL: Hydraulic 180kg or electric 300kg

Castor Size: Quality 125mm twin wheel castors with total lock and steer

Options: Both junior and wider sized available

Product Code

MST-SA-CTU

