

Contour Endo-X

The endoscopy imaging stretcher with a 'C' base for complete freedom of camera movement.



Key Benefits

Optimal design: Increased column spacing with left or right hand chassis

Supercharged: Keep things moving with long-life, cordless, fast-charge battery packs that work all day.

Custom quality: Manufactured in Australia with speed, quality, and care, plus customisation capability.

Longest warranty: Peace of mind with a lifetime of local support, easy maintenance, and fast guaranteed parts.

Application

Ideal for Endoscopy & Imaging

Why Buy?

An increased distance between the lifting columns of this stretcher allows excellent coverage with the C-Arm X-Ray camera, avoiding the risks of re-positioning and moving patients during procedures. Complete freedom of movement is also facilitated by the single rear bar joining the stretcher base, which creates ample space for the camera and leaves no restriction to the working height required by staff.

v-08122023

Contour Endo-X



Base option



Ultra-low height

Very safe and easy for patients to get on and off.

The 570mm low platform height of this table takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Excellent high height

Preserve your back and raise a patient up to your best working height.

Lifting the platform height up to 940mm (plus an extra 100mm for the mattress pad thickness), helps you to examine and treat your patient without having to stoop over, ensuring easy access and less stress. The backrest travels from 0 to 85 degrees with emergency CPR quick release.



Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



Easy customisation

Increase top width, change colours, mattress and dropside styles.

Whether it is a specialised medical clinic or a department in a large hospital, every single customer has their own unique requirements to be able to function safely and efficiently. This stretcher has a longer, 7 riser PVC upgrade dropside.

Specifications

Overall Size: 2135 x 815mm (705mm narrow & 965mm wide options)

Patient Surface: 1900 x 680mm (570mm narrow & 830mm wide options)

Height Range: 570 - 940mm

SWL: 300kg

Castor Size: 200mm large bearing castors

Package Includes: 100mm premia mattress, folding IV, dropdown footboard, vertical oxy, split drop sides, emergency CPR backrest release, C-shaped base for maximum camera clearance, 2nd handset & USB port

Options: Split dropside option (shown) and mid-position catch upgrades

Product Codes

Contour Endo-X: MPT-C-EX

Contour Endo-X Package: MPT-C-EX-P

Split fold dropside option: MPT-CA-SDS

v-08122023