



Eye Hospitals



Eye Hospitals


It can be challenging to balance surgeon requirements, carer functionality and patient comfort. Whatever the procedure, our products are carefully designed to maintain this balance, allowing necessary patient access and effective use of space.



Contour Recline E-Vertex



Foot Board



Meal Tray



Arm Board



Padded Armrests



Color Options



Contour E-Vertex



Stretcher IV Pole



Oxy Holders



Dropdown Foot Board



Document Holder

Contour Recline E-Vertex

A world-first, fully electric procedure chair with touch control adjustable headrest for ease of operation, access and patient comfort.



Key Benefits

Full flexibility: The infinitely adjustable headrest, range of chair heights with lie-flat and forward/reverse options places the surgeon and patient in exactly the right position

Safety and support: Increased leg clearance, with an adjustable, rigid locking headrest and wrist support accessory, brings patients close to protect surgeons from neck and back strain.

Comfort and control: Instantly increase patient comfort and enhance access with one-touch, electronic control of the scalloped or flat cushioned headrest; and pillow support accessory

Hygiene conscious: Smooth, gap free platform aids cleaning and helps prevent fluid ingress to meet the highest infection control standards.

Application

Ideal for Eye, ENT, plastics and dental care.

Why Buy?

The Contour E-Vertex medical procedure chair has been designed with direct feedback and insights from the country's leading Eye, ENT, plastics and dental professionals. The breakthrough design combines state-of-the-art technology and ergonomics to deliver exceptional functionality and support for procedures requiring close proximity, control for the practitioner and patient comfort. With convenient, touch-button handsets and user-friendly mechanisms, allowing a surgeon and their support team to quickly, safely and smoothly adjust the headrest, chair height and position for the procedure, patient transport and recovery.

[Also available as a Stretcher, view the Contour E-Vertex here.](#) *International patents pending

Contour Recline E-Vertex



Low and safe

Excellent 20 inch low height with fold-under drop-sides.

The ultra-low height is combined with drop-sides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patient's legs or clothes.



Lies flat

The backrest and leg rest move independently or together from chair to bed.


You now have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat.



Smart transfer

Conveniently stored handles unfold for staff friendly transportation.

Lock in the double steer system, unfold the push handles and move forward effortlessly on the large castors.




Slim and clean

The ultra thin design safely hides all mechanical operation and has a one-piece, gap free top surface for easy cleaning.




Effortless control

A clearly marked handset infinitely positions the headrest with the touch of a button.



Raise and lower

Adjust the height up and down while remaining level.



High to low tilt

Excellent tilt adjustment at any height.



Wrist Support option

Coated, 11 inches wide, height adjustable and removable rest

Specifications

Size: 30.7 x 45.3 inch overall (68.9 inch lying flat)

Seat width: 25.6 inch wide seat with 4 inch wider option available and 4 inch premium mattress

Height Range: 20 - 34 inch height range with lie-flat and forward / reverse tilt

SWL: 660lbs

Castor Size: Twin wheel 6 inch castors with double steer

Package Option: Standard IV, vertical oxy holder, drop down footrest, meal tray, standard scalloped or flat headrest cushion option

Product Codes

Contour Recline Vertex; MPC-C-REVE

Wrist Support; MPT-CA-OWS

Pillow Support; MPT-CA-CO-OSS

Flat Headrest Cushion Option; E0311

Foot Board

A removable board that locks into a padded foot rest or folds out for extra length as a leg support.

Key Benefits


- Covered in heavy duty vinyl
- Fits standard Recline leg support
- Use as a foot rest or fold down to lengthen the leg support

Product Code

- MPC-CA-DDF

Why Buy

Clip in this foot board to provide a comfortable padded rest for the patients feet to help restrain them whilst moving the chair. It is covered in hard wearing vinyl that is easy to clean. The board can be folded down if the legrest is raised to provide extra length and support for the patient legs.



Meal Tray

Tough tray constructed from medical grade plastic, molded with cup and meal positions.

Key Benefits


- Can be fitted to all four corners of the chair
- Swivels away for ease of access
- Folds down for convenient storage

Product Code

- MPC-CA-MLT

Why Buy

Quickly add this compact meal tray onto either side of the chair, lift and lock into position and swivel around to suit the patient. The tray folds down for easy storage and is constructed of a tough, one piece medical grade plastic that is very easy to clean. There are recesses to help prevent items from sliding off. This tray can be customised to any shape or size if required.



Arm Board

Lightweight armboard that clips onto the dropside hand rail and is padded for patients comfort.

Key Benefits

- 25mm high density foam covered in PU material
- Solid medical grade plastic base board
- Clips on both inside and outside of the dropside
- Adjustable angles for best patient positioning

Product Code

- MPC-CA-SAB

Why Buy

Simply and securely lock the arm board in a choice of locations ranging from the front to the back of the dropside that suits the size and position of the patient. It may also be locked into a variety of angles to help the carer's access. The design is light to carry and store, quick to use and easy to clean.



Padded Armrests

Thick clip-on pads to soften the dropside armrests.

Key Benefits

- High density foam core with slow recovery outer
- Covered in welded, pressure relieving stretch material
- 14.5 long x 4.9 wide x 8.3 inches high
- 1.6 inches of padding on top

Product Code

- MPC-CA-PDA

Why Buy

The Contour Recline Chair has been designed as a practical, clinical product that is functional and easy to clean. In some cases a patient may require extra support or comfort if spending longer periods of time in the chair. These easy to handle and store pads clip on over the dropside to widen and soften the armrest surfaces.



Color Options

A selection of colors to customize the Contour Recline procedure chairs for easy identification.



Key Benefits

- Tough extruded covers for a comfortable hand grip and bump protection.
- Used on chair dropside hand rails.
- Easy to replace for colour change or wear and tear.

Color Options



Flame Red



Midnight Blue



Deep Violet (Purple)

Product Codes

- Flame Red - MPT-CC-FR
- Midnight Blue - MPT-CC-MB
- Deep Violet (Purple) - MPT-CC-DV

Contour E-Vertex

A fully electric stretcher with effortless head positioning and comfortable access.



Key Benefits

Full flexibility: The infinitely adjustable headrest up and down tilt and excellent raise height allows the surgeon to place a patient in exactly the right position.

Safety and support: The compact, rigid locking headrest with increased leg clearance brings patients close and keeps them stable while minimizing surgeon neck and back strain.

Comfort and control: Instantly increase patient comfort and enhance access with one-touch, electronic control of the scalloped or flat cushioned headrest; and lightweight removable frame with cushions to support a pillow during recovery.

Slim and clean: The ultra-thin design hides all mechanical operation and has a one-piece, gap-free top surface for ease of cleaning and infection control.

Application

Ideal for Eye, ENT, plastics and dental care.

Why Buy?

We know precision care requires 100% confidence in the equipment being used. That's why the Contour E-Vertex combines state-of-the-art technology and ergonomic design to deliver exceptional support for the performance of procedures requiring close patient proximity, control for the practitioner and patient comfort. The slimline design of this extremely stable, rigid locking headrest provides excellent leg clearance – bringing the patient closer. While the touch button electronic positioning system* allows a surgeon or their support team to quickly, safely and precisely position a patient's head. [Also available as a procedure chair, see the Contour Recline E-Vertex here.](#) *International patents pending

Contour E-Vertex



Slim and clean

The ultra thin design safely hides all mechanical operation and has a one-piece, gap free top surface for easy cleaning.



Effortless control

A clearly marked handset infinitely positions the headrest with the touch of a button.



Raise and lower

Adjust the height up and down while remaining level.




High tilt

Excellent raise height with up and down tilt adjustment.



Low tilt

A practical low height with up and down tilt adjustment.



Pillow Support option

Lightweight removable frame with cushions to help support a pillow during recovery.



Wrist Support option

Coated, 11 inches wide, height adjustable and removable rest.



Fold-back design

The wrist support flips away to allow quick and easy access to the patient.

Specifications

Size: 87 x 32 inches overall with 38 inch wide option

Height Range: 20 - 34 inches

SWL: 660Lb

Castor Size: 8 inch castors and Ezi-Track 5th wheel steering

Package Includes: 75mm premia mattress, vertical oxy holder, folding IV, USB Port, 2nd Handset, Emergency CPR Release

Product Codes

Contour E-Vertex; MPT-C-VE-E

Wrist Support; MPT-CA-OWS

Pillow Support; MPT-CA-CO-OSS

Flat Headrest Cushion; MPT-CA-FHC

Stretcher IV Pole

Include a standard removable or fold down IV Pole with your stretcher.

Key Benefits

- A strong and safe hook design
- Light and easy 'lift out' auto-lock hanger
- Non-rotational design to prevent tube tangling
- High strength locking for heavy loads
- Tough construction from stainless steel and nylon
- Easy 'pull-down' release from any angle
- Large 1.1 inch diameter stand body for ease of pump mounting
- Friction retainer mount for effortless use
- Removable and permanent fold down options
- 28 inch up to 44 inch total height
- 22Lb weight capacity, 11Lb per hook
- Only 2.2Lb total weight

Product Codes

- Standard IV Pole - MPT-CA-IVP

Why Buy

Constructed tough from stainless steel and nylon, these clever non-rotational IV poles are designed to withstand the harshest use and make daily life easier by helping to prevent the tangling of critical soft tubing and cords. Simply lift out the hanger to any height and it automatically locks into an infinite position. Pull the lock body down from any direction to release and lower. The standard removable IV fits into all four corners of a stretcher top and stores underneath in the base if there is provision. The folding IV fits onto either end of the stretcher top and can be installed to fold down either to the left or right.



Standard IV storage



Folding IV Pole

Oxy Holders

A selection of stainless steel oxy bottle holders that can be mounted vertically or horizontally.

Key Benefits

- Vertical basket size for up to 4.7 inch diameter bottle; 12.6 deep x 6 inches wide
- Horizontal basket size for up to 6.3 inch diameter bottle; 23.6 long x 7 inches wide
- Horizontal basket can vary in hanging height depending on stretcher type; typically 11 inches

Product Codes

- Vertical Rail Mount - MPT-CA-VO2
- Vertical Corner Mount - MPT-CA-VCO2
- Vertical Corner Upright Mount - MPT-CA-VCUO2
- Horizontal Oxy Holder - MPT-CA-HO2

Why Buy

Due to so many wide and varied departmental requirements for specific oxy bottle sizes and positioning, we have developed a range of options. Larger bottles can be laid horizontally to ensure they are kept out of the way and smaller bottles are commonly mounted vertically. Some prefer 'fixed', non-moving positions and others require 'swinging' options that allow the bottles to be moved aside if they become an obstruction during transport. There is also the option to support bottles vertically either above or below the patient surface.



Vertical Corner Mount



Vertical Corner High Mount



Horizontal Oxy Holder



Horizontal Oxy Holder

Dropdown Foot Board

Padded fold-down foot board to increase the length of the patient surface.

Key Benefits

- High density foam covered in welded stretch material
- Tough board construction in a one-piece hospital grade plastic
- Stainless steel lift and lock folding hinge
- 21 wide x 13 inches high (10 inch high cushion)
- Adds an extra 12 inches to the length of a mattress

Product Code

- MPT-CA-DDF

Why Buy

Padded for patient comfort, the removable dropdown foot board has molded handles for pushing the stretcher and with a simple 'lift and lay' action, folds down to quickly increase the length of the mattress area.



Fold down

Document Holder

Strong powder coated wire basket that is hung or mounted for document storage.

Key Benefits

- Ultra-tough, powder coated, 0.19 of an inch steel wire construction
- External size; 13 inch wide x 3 inch deep x 9 inch high
- Internal size; 15 inch wide x 2 inch deep x 9 inch high
- Hang on the inside or outside of a stretcher end rail
- Wall and accessory mounting options

Product Code

- MPT-CA-DCH

Why Buy

Document holders can be used in a variety of ways for many different storage requirements. In the past we have mounted holders on either end or in special places around a stretcher, off accessories such as the dropdown foot board and even on the walls around the stretcher holding bays. We can customize sizes and designs and will work with you to ensure you can maximize department efficiency.



Removable holder

Notes

Notes
