

Dropdown Foot Board

Padded fold-down foot board to increase the length of the patient surface.

Key Benefits

- High density foam covered in welded stretch material
- Tough board construction in a one-piece hospital grade plastic
- Stainless steel lift and lock folding hinge
- 21 wide x 13 inches high (10 inch high cushion)
- Adds an extra 12 inches to the length of a mattress

Product Code

- MPT-CA-DDF

Why Buy

Padded for patient comfort, the removable dropdown foot board has molded handles for pushing the stretcher and with a simple 'lift and lay' action, folds down to quickly increase the length of the mattress area.



Fold down

v-23012023